



**STARKSTROM**

disposable  
scope warmer

Quality Medical Equipment and Service



DOVIC OUT-2026, YE  
230  
PELUPS 2020N ALASH TO A  
UPD 2020N ALASH 147



# scope warmer

when clarity of vision  
really matters



## disposable scope warmer

Eliminates Fogging During Laparoscopic Procedures



### Quick & easy to use - Sterile packaged Heats up the scope(s) in one minute

- | No batteries; no electricity
- | No capital equipment required
- | No alcohol
- | Environmentally safe
- | CE marked
- | Scope stays warm for 45 minutes to one hour
- | Backup scope is always heated
- | Accommodates all sizes and lengths of scopes
- | Simple one-step activation
- | Self-contained instant heat source



256 Field End Road Eastcote Ruislip Middlesex HA4 9UW

T: +44 (0) 208 868 3732 F: +44 (0) 208 868 3736

E: [info@starkstrom.com](mailto:info@starkstrom.com) [www.starkstrom.com](http://www.starkstrom.com)